

REVEALCOVID-19™



Total Antibody Test

medmira
helping people know...®

REVEALCOVID-19® PLUS

REVEALCOVID-19® PLUS is a manually performed, visually interpreted rapid test designed for the qualitative detection of antibodies against the nucleocapsid (N) and spike (S) proteins of SARS-CoV-2, the causative agent in the COVID-19 disease, allowing the detection of the current variants of interest (Delta, Alpha, Beta).

About MedMira

MedMira's patented Rapid Vertical Flow™ Technology delivers fast, accurate results in a wide range of applications, it only takes few minutes and three simple steps to go from sample to result. Our products help healthcare providers and their patients know more, instantly.



Technology & Performance

REVEALCOVID-19® PLUS does not require specialized equipment or trained laboratory personnel to be performed. The membrane is coated with recombinant antigens corresponding to the regions of SARS-CoV-2 nucleocapsid and spike proteins.

When used to test SARS-CoV-2 positive samples, REVEAL COVID-19® PLUS has a positive percent agreement (PPA) of 99.02% (n=102; 8+ days post symptom onset) and a negative percent agreement (NPA) of 97.78% (n=451) using serum and plasma samples. The combined specificity of the test was calculated to be 100% (n=89).

Features

- Simple test procedure
- Instant and easy-to-interpret results
- Works with whole blood, serum, and plasma
- No timers or specialized equipment required
- High sensitivity and specificity
- Current variants (Delta, Alpha, Beta) can be detected

Ordering Information

POU (Fingerstick Whole Blood) [Cat. No. 815211006357; Quantity/Box: 20]

Applications

- Pharmacies and Physician clinics
- Occupational testing sites
- Prevalence testing
- Public health strategy planning